

## **NIDB Medium 13**

### **Potato Dextrose Agar (PDA)**

#### **Agar Medium:**

Potato Dextrose Agar (BD213400) ..... 39.0 g  
DI Water ..... 1000 mL

Mix thoroughly

Autoclave at 121°C for 15 mins.

\* This product is commercially available from BD.

To make medium from scratch, follow formulation below:

Potato Starch (from infusion) ..... 4.0 g  
Dextrose ..... 20.0 g  
Agar ..... 15.0 g  
DI Water ..... 1000 mL

Final pH 5.6 ± 0.2

Mix thoroughly

Autoclave at 121°C for 15 mins.