

NIDB Medium 32

Yeast extract-peptone-dextrose (YPD) Medium (Broth)

Broth Medium:

YPD Broth (BD 242820) 50.0 g
DI Water 1000 mL

Mix thoroughly

Autoclave at 121°C for 15 mins.

* YPD is commercially available from BD.

To make medium from scratch, follow formulation below:

Yeast Extract 10.0 g
Peptone 20.0 g
Dextrose 20.0 g
DI Water 1000 mL

Final pH 6.5 ± 0.2

Mix thoroughly

Autoclave at 121°C for 15 mins.