

NIDB Medium 33

Yeast extract-peptone-dextrose (YPD) Medium (Agar)

Agar Medium:

YPD Agar (BD 242720)	65.0 g
DI Water	1000 mL

Mix thoroughly

Autoclave at 121°C for 15 mins.

* YPD is commercially available from BD.

To make medium from scratch, follow formulation below:

Yeast Extract	10.0 g
Peptone	20.0 g
Dextrose	20.0 g
Agar	20.0 g
DI Water	1000 mL

Final pH 6.5 ± 0.2

Mix thoroughly

Autoclave at 121°C for 15 mins.